

Analysis of Saccharides/ 糖類の分析

The NH₂ column is often used for analysis of saccharides, and separation shows similar behavior to HILIC. We analyzed saccharides using the Develosil ANIDIUS, a column that is very durable even in highly mobile phase.

糖類の分析にはNH₂カラムを用いることが多く、分離はHILICと同じ挙動を示します。そこで水の多い移動相においても耐久性のある「Develosil ANIDIUS」を用いた糖類の分析例を紹介します。




Fig.1 Analysis of Saccharides

Conditions;

Column: Develosil ANIDIUS
(4.6x250mm)

Mobile phase: Acetonitrile/Water=80/20

Flow rate: 1.0ml/min

Temperature: 40°C

Detection: RI

Sample: Fructose, Glucose, Sucrose,
Maltose

Injection volume: 20uL

System: 日本分光製 / JASCO 200




Fig.2 Analysis of Maltooligosaccharides

Conditions;

Column: Develosil ANIDIUS
(4.6x250mm)

Mobile phase: Acetonitrile/Water=65/35

Flow rate: 1.0ml/min

Temperature: 40°C

Detection: RI

Sample: Glucose, Maltose, Maltotriose,
Maltotetraose, Maltopentaose,
maltohexaose, Maltoheptaose

Injection volume: 20uL

System: 日本分光製 / JASCO2000