


S. Horikiri, Nomura Chemical | 10060 Carroll Canyon Rd, Ste 100, San Diego CA 92131

■ Introduction

High resolution chromatography is a necessity to ensure consistent and accurate results. In order to fulfill this need, FlexFire is introducing the FlexFire Wide Pore (WP) C18 1.6 µm (2.0 x 50 nm) column. The new 1.6 µm particle size offers higher resolution, compared to 2.6 µm particle sizes, while maintaining the great resolution for a 1 base difference in RNA.


In addition to analysing and displaying a side by side comparison between the FlexFire WP C18 1.6 µm (2.0 x 50 nm) and the FlexFire WP C18 2.6 µm (2.0 x 50 nm), Separation analysis comparing the mobile phase additives of DIPA to EDIPA have been provided.

■ Comparison by particle size with FlexFire WP C18 2.6 µm (2.0 x 50 nm)




Analytical Conditions					
Column	FlexFire WP C18, 1.6 µm (2.0 x 50 mm)				
	FlexFire WP C18, 2.6 µm (2.0 x 50 mm)				
Mobile Phase	A) 100 mM HFIP + 10 mM DIPA				
	B) A/Methanol = 50/50				
Gradient 1.6 µm	min	mL/min	%A	%B	Curve
	0.00	0.3	80	20	
	5.04	0.3	40	60	6
	5.05	0.3	80	20	6
Gradient 2.6 µm	min	mL/min	%A	%B	Curve
	0.00	0.3	80	20	
	8.40	0.3	40	60	6
	8.42	0.3	80	20	6
Temperature	60 °C				
Detection	UV 260 nm				
Sample	40-mer miRNA				
Injection Volume	1.0 µL				
System	Waters ACQUITY UPLC H-Class PLUS				

Mobile Phase Additives With WP C18 1.6 μm



Analytical Conditions					
Column	FlexFire WP C18, 1.6 μm (2.0 x 50 mm)				
Mobile Phase	A) 100 mM HFIP + 10 mM DIPA				
	B) A/Methanol = 50/50				
Mobile Phase	A) 100 mM HFIP + 10 mM EDIPA				
	B) A/Methanol = 50/50				
Gradient 1.6 μm	min	mL/min	%A	%B	Curve
	0.00	0.3	80	20	
	5.04	0.3	40	60	6
	5.05	0.3	80	20	6
Temperature	60 ° C				
Detection	UV 260 nm				
Sample	40-mer miRNA				
Injection Volume	1.0 μL				
System	Waters ACQUITY UPLC H-Class PLUS				

Resolution For a 1 base Difference in RNA



Analytical Conditions					
Column	FlexFire WP C18, 1.6 μm (2.0 x 50 mm)				
Mobile Phase	A) 100 mM HFIP + 10 mM DIPA				
	B) A/Methanol = 50/50				
Mobile Phase	A) 100 mM HFIP + 10 mM EDIPA				
	B) A/Methanol = 50/50				
Gradient 1.6 μm	min	mL/min	%A	%B	Curve
	0.00	0.3	80	20	
	5.04	0.3	40	60	6
	5.05	0.3	80	20	6
Temperature	60 ° C				
Detection	UV 260 nm				
Sample	1: GUACGCGGAUACUUCGAtt (20mer) 2: CGUACGCGCAUACUUCGAtt (21mer)				
Injection Volume	1.0 μL				
System	Waters ACQUITY UPLC H-Class PLUS				

■ Contact Us



Develosil USA

10060 Carroll Canyon Rd. Ste. 100
San Diego, CA 92131
Phone: 858-800-2433
Web: <https://develosil.us/>



Nomura Chemical Co., Ltd.

15, Hinode-cho, Set, 489-0004, Japan
Tel: +81-561-48-1853
Fax: +81-561-48-1434
E-mail: info@develosil.net