




■ Calibration curve of 40 mer miRNA




The calibration curve obtained by the UV detector presented a clean straight line in the range of 0.32-650 ng of miRNA. LC/MS should yield similar results at even lower concentrations.

■ Detection of 1-base difference in miRNA



Resolution has been improved in the FlexFire wide pore columns, and the FlexFire Wide Pore (WP) C8 has the same retention as the FlexFire WP C18. The performance of these columns is dramatically improved compared to the conventional product.



Analytical Conditions

Column

- FlexFire C18, 2.6 μm (2.0 x 50 mm)
- FlexFire WP C18, 2.6 μm (2.0 x 50 mm) - Metal free
- FlexFire WP C18, 2.6 μm (2.0 x 50 mm) - Stainless
- FlexFire WP C8, 2.6 μm (2.0 x 50 mm)
- FlexFire WP C4, 2.6 μm (2.0 x 50 mm)
- FlexFire mAb-RP, 2.6 μm (2.0 x 50 mm)

Mobile Phase

- A) 100 mM HFIP + 10 mM TEA
- B) Mobile Phase A/Methanol=50/50

Gradient	min	mL/min	%A	%B	Curve
	0.00	0.3	93	7	
	8.60	0.3	73	27	6
	8.61	0.3	93	7	6

Detection UV 260 nm


Sample 20 nt (61.4 $\mu\text{g}/\text{mL}$)
21 nt (63.9 $\mu\text{g}/\text{mL}$)

Injection Volume 2.0 μL

System Waters ACQUITY UPLC H-Class PLUS

Mixer Volume 100 μL

■ Effect of HFIP concentration



Analytical Conditions

Column	FlexFire WP C18, 2.6 μ m (2.0 x 50 mm)				
Mobile Phase	A) 200 mM HFIP + 20 mM TEA				
	A) 100 mM HFIP + 10 mM TEA				
	A) 50 mM HFIP + 5 mM TEA				
	B) Mobile Phase A/MeOH+50/50				
Gradient	min	mL/min	%A	%B	Curve
	0.00	0.3	93	7	
	8.60	0.3	73	27	6
	8.61	0.3	93	7	6
Temperature	60°C				
Sample	1. 5'-GUACGCGGAAUACUUCGAdTdT-3' (20 nt)				
	2. 5'-CGUACGCGGAAUACUUCGAdTdT-3" (21 nt)				
Injection Volume	1.0 μ L				
System	Waters ACQUITY UPLC H-Class PLUS				
Mixer Volume	100 μ L				

HFIP concentration and retention.

HFIP concentration was altered while the ratio of [HFIP] : [Ion Pair (base)] was kept at 10:1. At 5 mM HFIP, the peak appeared in the front peak. At 50 mM HFIP, a sharp peak was detected. However, the peaks of 20 nt and 21 nt tended to overlap slightly.

Increasing HFIP concentration to 100 mM showed good separation. At 200 mM, the separation power was better but the effect of HFIP seemed to have plateaued. We recommend using 100 mM HFIP as it shows good separation performance while minimizing damage to the column.

■ Conclusion

Important points when selecting a column for the analysis of nucleic acid oligomers are:

- Using a wide pore column
- Use a silica gel substrate with very low impurities

In addition to HFIP, TEAA (triethylammonium acetate) and ion pair concentrations affect the separation. When using LC/MS, the mobile phase constituents can affect both the sensitivity and pollution of the detector.

■ Contact Us



Develosil USA

10060 Carroll Canyon Rd. Ste. 100
San Diego, CA 92131
Phone: 858-800-2433
Web: <https://develosil.us/>



Nomura Chemical Co., Ltd.

15, Hinode-cho, Set, 489-0004, Japan
Tel: +81-561-48-1853
Fax: +81-561-48-1434
E-mail: info@develosil.net