


Analysis of bromhexine hydrochloride/ブロムヘキシン塩酸塩の分析



Sample

1.o-Hydroxybenzamide (0.15mg/ml)

2.Bromhexine Hydrochloride (0.15mg/ml)

Conditions;

Column: Develosil XG-C18M (5um) 4.6x150mm

Mobile phase: Acetonitrile/50mM Ammonium Formate=80/20

Flow rate: 1.0ml/min

Temperature: 40°C

Detection: UV254nm

Injection volume: 2.0uL

System: DIONEX/Ultimate 3000 Series