


局方：エフェドリン塩酸塩の分析/Analysis of Ephedrine hydrochloride



Conditions;

Column: Develosil XG-C18LC-5

Size: 4.6x250mm

Mobile phase: ラウリル硫酸Na(1→128)/アセトニトリル/リン酸=640/360/1
Na lauryl sulfate (1 → 128) / acetonitrile / phosphoric acid = 640/360/1

Temperature: 45°C

Detection: UV210nm

System: 日本分光/LC-2000plus series

Injection volume: 10uL