


## Develosil Application Data

### No.071 Ingredients in cold medicines/かぜ薬成分



#### Analytical Conditions:

Column: Develosil HSR C18,5um (4.6x150mm)

Mobile phase: Methanol/50mM NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub>, pH2.5=40/60

Flow rate: 1.0mL/min

Temperature: 40°C

Detection: UV280nm

Sample: 1.Acetaminophen (0.13mg/mL)

2.Caffeine (0.03mg/mL)

3.Salicylamide (0.14mg/mL)

4.Aspirin (0.58mg/mL)

5.Ethenzamide (0.28mg/mL)

Injection volume: 5.0uL

System: Jasco 2000series