


Develosil Application Data

No.060 Local anesthesia/局所麻醉剂



Analytical Conditions:

Column: Develosil HSR C18,5um (4.6x150mm)

Mobile phase: Acetonitrile/25mM (NH₄)₂HPO₄, pH7.0=50/50

Flow rate: 1.0mL/min

Temperature: 40°C

Detection: UV210nm

Sample: 1.Procaineamide (0.11mg/mL)

2.Mepivacaine HCl (0.11mg/mL)

3.Tetracaine HCl (0.10mg/mL)

4.Lidocaine (0.12mg/mL)

5.Bupivacaine HCl (0.09mg/mL)

Injection volume: 2.0uL

System: Waters alliance